

Important safety information



Warning!

Follow all health and safety regulations.

- Use appropriate personal protective equipment
- Handle all samples as if they contain infectious agents.
- Observe all precautions and warnings in the Instructions for Use.

Document information

Document version: Initial version

L29SC3REN0

Issue date: 2020.08

Revision date: September, 2020



SD Biosensor, Inc.
C-4th&5th, 16,

Deogyong-daero 1556beon-gil,

Yeongil-gu, Suwon-si,

Gyeonggi-do, 16690,

REPUBLIC OF KOREA

Made in Korea

Distributed by:

Roche Diagnostics GmbH

Sandhofer Strasse 116

D-68305 Mannheim, Germany

www.roche.com

For in vitro diagnostics use only
Not for self testing

© 2020. All rights reserved.

IVD



SARS-CoV-2 Rapid Antigen Test Quick Reference Guide


This guide is a reference for using the
SARS-CoV-2 Rapid Antigen Test.
Read the Instructions for Use before using this test.

0 9327762001 (01) 2020-09 EN

1 Preparing for a test

1. Carefully read the Instructions for Use for the
SARS-CoV-2 Rapid Antigen Test.



2. Check the expiry date  on the back of the foil pouch. Do not use the test if the expiry date has passed.



3. Open the foil pouch and remove the test device and the desiccant package.

4. Ensure that the test device is undamaged and that the desiccant status indicator shows valid (yellow).



5. Perform a QC as required according to the Instructions for Use of the QC material and according to your local guidelines.

2a Collecting and preparing a sample (Nasopharyngeal swab)

1. Insert a sterile swab into the patient's nostril, swab

2. Insert the swab into an extraction buffer tube.

3. Remove the swab while squeezing the sides of the

4. Press the nozzle cap tightly onto the tube.